



CO₂RE[®] Intima
Women's Intimate Wellness



Women's
intimate wellness



Relieves
signs of childbirth & aging



In-demand
new treatment

Women's Intimate Wellness – An In-Demand, Rapidly Growing Need

Women's intimate wellness is a rapidly growing need.

In a survey of American women¹:

- 93% had not had any vaginal rejuvenation treatments.
- 43% had heard of **vaginal rejuvenation**, labiaplasty, labial tightening and pigmentation treatment.
- 63% would consider a non-surgical procedure to **tighten the vaginal canal**.
- 75% suffer from **vaginal dryness** and of these, 90% say that it **negatively affects their lives**.
- 66% suffer from **Stress Urinary Incontinence** to some degree.

The women's intimate wellness market is growing and women are taking control.



Many women experience painful intercourse, dryness, itching and burning of the vagina.²

Women have a growing interest in modifying the appearance of the vulva, including pigmentation.

With CO₂RE Intima, physicians are now able to offer treatment of these conditions.

CO₂ lasers have been shown to remodel tissue with neof ormation of collagen and elastic fibers in atrophic skin.³



Relieves signs of childbirth & aging

- In a survey of American women, 63% felt **their vagina was looser** since giving birth.
- Of these women, 62% felt **this negatively affected their sex life**.



In-demand new treatment

The emerging area of women's intimate wellness creates new treatment opportunities for health care providers. Educational workshops are available to guide you through skills acquisition.



Women Report High Satisfaction with CO₂RE Intima Treatments⁴

- 81% reported an improvement in **sexual gratification**.
- 94% reported an improvement in **vaginal rejuvenation**.
- 100% reported **satisfaction with treatment** and 94% would recommend it.

CO₂RE Intima Uniquely Treats All 3 Intimate Areas: Internal, External and Introitus.



Treat vulvovaginal atrophy (VVA) and do vaginal rejuvenation. Remodel vaginal connective tissue³ with the hygienic, single-use Internal handpiece.



Treat the introitus and improve the appearance of the external genitalia with the hygienic, single-use External handpiece.



Use the Surgical handpiece to excise or incise external tissue.



Women's
intimate wellness

CO₂RE Intima provides a value packed solution, attractively priced for easy access to this emerging women's intimate wellness market.

With CO₂RE Intima, the hygienic, single-use handpieces assure women of the hygiene they expect of treatments in this intimate area.

These treatments are easy-to-master and quick with both internal and external treatments taking 10 minutes.

non-surgical  in-office  little to no downtime  mild to no discomfort

CO₂RE Intima leverages the CO₂RE platform

CO₂RE Intima is the latest enhancement to the versatile CO₂RE platform which offers proven aesthetic solutions for skin resurfacing, wrinkle and scar reduction, minor dermal incisions, and now women's intimate wellness treatments.

CO₂RE gives you the ability to treat both superficial and deep skin layers simultaneously with precision control over the intensity, pattern and depth of ablation – all from a single platform, letting you differentiate and efficiently grow your practice.

- 4 Modes: Light, Mid, Deep & Fusion for fractional resurfacing
- Classic Mode for full resurfacing
- Surgical Mode for incision

- Single-use Internal & External handpieces and lenses
- Multi-use Surgical handpiece and lens
- Easy-to-learn user interface
- Lightweight articulated arm for ease-of-use
- Compact size & easily mobile



CO₂RE System Specifications

| | |
|---------------------------|--|
| Laser Type | CO ₂ |
| Wavelength | 10,600 nm |
| Laser Emission Mode | Pulsed, CW |
| Output Power | 60 Watts (maximum peak) |
| Laser Beam Energy | Up to 80 mJ |
| Pulse Repetition Rate | Up to 16.7 kHz |
| Pulse Duration | Up to 1100 µsec |
| Maximum Scan Area | 10 mm diameter |
| Aiming Laser | 5 mW, 650 nm (red) diode laser |
| Articulated Arm | <ul style="list-style-type: none"> • Lightweight, aluminum fiber, 7-joint, pneumatically balanced • Treatment radius at full arm extension: 85 cm (33.5") • Horizontal rotation: 360° |
| Lens Assemblies | Two 150 µm spot size (1 swivel). 120 µm spot size. |
| Handpieces | Hygienic, single use Internal and External handpieces. Surgical handpiece is multi-use. |
| User Interface | 10.4" high-resolution color LCD Touch-screen technology |
| Laser Emission Control | Footswitch |
| Laser Emission Indicators | Solid amber LED in Standby and Ready modes. Flashing red LED during laser emission. |
| Electrical Requirements | 100 - 240 VAC self-sensing, 50/60 Hz, 7.1 A, single phase. |
| Dimensions | 40 cm x 44 cm x 177 cm (15.75" x 17.3" x 70") |
| Weight | 26 kg (58 lbs) |

Operating Modes

| Mode | Energy Settings (mJ) | Fluence (J/cm ²) | Fractional Coverage (%) | Estimated Ablation Depth (µm) |
|-----------------------------|----------------------|------------------------------|-------------------------|---|
| CO ₂ RE Light | 18 - 74.5 | 3.6 - 6.6 | 30 - 50 | 20 - 40 |
| CO ₂ RE Mid | 45 - 257 | 8.8 - 22.7 | 20 - 40 | 100 - 200 |
| CO ₂ RE Deep | 30 - 80 | 170 - 453 | 1 - 5 | 500 - 750 |
| CO ₂ RE Fusion | 45 - 128 30 - 80 | 8.1 - 24.8 170 - 453 | 20 - 40 | 100 - 200 & 500 - 750 |
| Classic Resurfacing | 1 - 10 | 5.6 - 56.5 | Full | 100 - 200 |
| Dermal Excision or Incision | 5 - 40 | 28.3 - 226 | N/A | 0.2 mm x selected line 0.2 mm freehand |

Women's Intimate Wellness - A Rapidly Growing Need



Science. Results. Trust.

Syneron Candela is a leading global aesthetic device company with a comprehensive product portfolio and a global distribution footprint.

The Company's technology enables physicians to provide advanced solutions for a broad range of medical-aesthetic applications including body contouring, hair removal, wrinkle reduction, tattoo removal, improving the skin's appearance through the treatment of superficial benign vascular and pigmented lesions, and the treatment of acne, leg veins and cellulite.

The Company sells its products under three distinct brands, Syneron, Candela and CoolTouch, and has a wide portfolio of trusted, leading products, including UltraShape, VelaShape, GentleLase, VBeam Perfecta, PicoWay, CO₂RE, CO₂RE Intima*, Profound and elōs Plus.

Founded in 2000, the company markets, services and supports its products in 86 countries. It has offices in North America, France, Germany, Italy, Portugal, Spain, UK, Australia, China, Japan, and Hong Kong, and distributors worldwide.

SYNERON  CANDELA®

www.syneron-candela.com

¹ Proprietary Syneron Consumer Survey, 2015

² Perino A. Vulvo-vaginal atrophy: A new treatment modality using thermo-ablative fractional CO₂ laser. *Maturitas*. 2015 Mar;80(3):296-301.

³ Salvatore S, et al. Histological study on the effects of microablative fractional CO₂ laser on atrophic vaginal tissue: an ex vivo study. *Menopause*. 2015, Jan 20. [Epub ahead of print].

⁴ Clinical Study HM Montepincipe Hospital, Madrid, Spain, 2015

* CO₂RE Intima is a module of CO₂RE